







8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • WWW.HOBBSMEDICAL.COM



HOBBS MEDICAL, INC.



INTRODUCTION

WELCOME

Established in 1984, Hobbs Medical, Inc.

has built a solid reputation for

innovation and quality. A privately

held company with expertise in the

development and manufacture of

diagnostic and therapeutic accessories for

flexible GI and Pulmonary Endoscopy, Hobbs

works with a network of industry professionals,

consultants, sales representatives and international distribution agencies

to bring precision endoscopic accessories to the world market. We take pride in

our comprehensive line of precision accessories, all manufactured in the USA, and

continue to introduce unique devices to the ever-changing market.

The entire staff at Hobbs Medical is committed to exceeding each and every customer's expectations;

providing cost-effective, quality and innovative products and investing in ongoing product development.

Our Sales personnel are always available to answer inquiries. We look forward to responding to your needs.

Contact Sales at (800) 344-6227 or visit www.hobbsmedical.com



HOBBS MEDICAL, INC. 8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • WWW.HOBBSMEDICAL.COM

ALL HOBBS MEDICAL PRODUCTS ARE PROUDLY MANUFACTURED IN THE USA.
HOBBS MEDICAL IS ISO-13485 REGISTERED; CE-CERTIFIED AND COMPLIES WITH
COUNCIL DIRECTIVE 93/94EEC AND FDA QUALITY SYSTEMS REGULATIONS.

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

QUALITY POLICY

The Quality Policy of Hobbs Medical, Inc. is to meet the needs and expectations of our customers. This will be accomplished through market feedback and continuous evaluation of our Quality Management System to ensure its continued effectiveness and suitability in complying with customer needs and all regulatory requirements at all levels of the organization.

SPECIFICATION CHANGES

Hobbs Medical reserves the right to make improvements and changes in the products described herein. Should you have a suggestion for a new device or an improvement of an existing device, let our Engineering staff translate your idea into reality. Our website (www.hobbsmedical.com) provides an avenue for product suggestions.

PRODUCT WARRANTY

Descriptions or specifications appearing in Hobbs Medical literature are solely meant to generally describe the product at the time of manufacture and do not constitute any express warranties. Hobbs Medical warrants all of its products to be free from defects in materials and workmanship. Storage conditions and product use are factors beyond the control of Hobbs Medical and may directly affect the products and the results obtained from them. The user assumes responsibility for any consequences, which may arise from the final use of the products. No person has any authority to bind Hobbs Medical to any representation or warranty concerning these products.

All reusable products in this catalog are warranted to the original purchaser to be free from defects in material and workmanship for a period of six months from the date of purchase. All disposable products are warranted to the original purchaser to be free from defects in materials and workmanship for the first use only.

RETURN POLICY

A Return Authorization Number must be obtained for all returned products. Contact our Customer Service Department (800-344-6227). Returns must include a reason for return, lot number, catalog number and invoice number. Contaminated products sent for evaluation will not be accepted unless the product has been cleaned and high level disinfected prior to returning to Hobbs Medical. Returned, unused products must be in the original unopened Hobbs Medical pouch and will be subject to a 20% restocking fee.

SHIPPING & TERMS

Most items are shipped within one day. Products are shipped via ground transportation unless otherwise specified. MasterCard and VISA accepted. F.O.B. Stafford Springs, Connecticut. Due to fluctuating operating costs, we reserve the right to change prices without prior notice.

ORDERING INFORMATION

Our Sales and Customer Service personnel are always available to answer inquiries and to meet our customers' specific needs. Order and Price Requests can be placed by telephone, fax or email.

CONTACT INFORMATION

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

Orders & Returns

Product Information
Pricing
Specification Changes




SPECIALTY ITEMS

SPECIALTY ITEMS

GASTROINTESTINAL PRODUCTS

Active Cords

- Reusable insulated cords
- Active Banana fitting (Male connector) fits all Hobbs Medical accessories
- Variety of fittings provided to meet your electrocautery machine requirements
- No minimum order requirements
- Packaged non-sterile



ACTIVE CORDS


Catalog Number	Description	Cautery Unit	Distal End
4520	Active Cord	Universal Fitting	Male
4521	Active Cord	Banana Fitting	Male
4522	Active Cord	Olympus Fitting	Male
4523	Active Cord	Cameron Fitting	Male



GASTROINTESTINAL PRODUCTS

The Mistifier: Spray Catheter: Disposable

- Used in chromoendoscopy for early detection of flat lesions and polyps
- Consistent 360-degree spray mist pattern
- Disposable to prevent cross-contamination concerns
- Compatible with power irrigators
- Packaged sterile for single use: 5 per box



THE MISTIFIER: SPRAY CATHETER: DISPOSABLE

Catalog Number	Description	Sheath	Length
2190	Spray Catheter	2.6mm	260cm

Aspiration Catheter: Disposable

- Teflon Catheter with universal Luer-lock fitting for ease of passage through working channel
- Available with permanent or removable Luer-lock fitting
- Spiral pattern of holes at distal tip
- Compatible with 2.8mm or larger channels
- Packaged sterile for single use: 5 per box



ASPIRATION CATHETER: DISPOSABLE

Catalog Number	Description	Sheath	Length	Luer-Lock
2181	Aspiration Catheter	2.3mm	240cm	Permanent
2182	Aspiration Catheter	2.3mm	240cm	Removable



GASTROINTESTINAL PRODUCTS

Scope Introducer

- Provides easy irrigation access through working channel
- Allows accessories to move smoothly through working channel
- Compatible with most currently manufactured endoscopes
- Packaged sterile for single use: 5 per box




SCOPE INTRODUCER

Catalog Number	Description	Outer Diameter	Inner Diameter
4650	Scope Introducer	3.07mm	2.48mm

GASTROINTESTINAL AND PULMONARY PRODUCTS

Vacu-lok Aspirating Syringe

- Universal Luer-lock fitting
- Use during Endoscopic Ultrasound/Fine Needle Aspirate for maximum negative pressure
- Use with inflation monitor* for inflation and evacuation of balloon dilators
- Ergonomic design provides a unique locking system
- Disposable to prevent cross-contamination
- Packaged sterile for single use: 5 per box
- Size 20ml



*See pages 6, 9, & 10 for Inflation Monitor


VACU-LOK ASPIRATING SYRINGE

Catalog Number	Description
3900	Vacu-lok Aspirating Syringe

GASTROINTESTINAL PRODUCTS

Flex-Ez Guide Wires

- Provides a cost-effective alternative when a generic Guide Wire is required to assist in placing Esophageal Stents, Over the Wire Balloon Dilators or Achalasia Balloons
- Teflon coating provides easy passage through the endoscope's working channel
- Packaged sterile for single use



FLEX-EZ GUIDE WIRES

Catalog Number	Description	O.D.	Style	Length
3406	Floppy Tip	.038in	Teflon Coated	260cm
3410	Floppy Tip	.035in	Teflon Coated	260cm
3412	Floppy Tip	.035in	Teflon Coated	400cm



CONTACT INFORMATION

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Orders & Returns

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

Product Information
Pricing
Specification Changes